



HOKKAIDO
UNIVERSITY

Pump-Probe **2 beam** XAS spectra

Hokkaido University
Kiyotaka Asakura

In order to understand Catalytic Reaction


Catalytic reaction.

To decide the active site and determine its structure.

To determine the intermediate structure and reaction step.

To determine the kinetics and the mechanism




To write the reaction cycle.




Wilkinson Catalyst

Correlation = simultaneous measurement



Ni2P



Proposal: Photoexcitation process and carrier




5\$5\$ 2\$ swlne\$NTGG\$459\$5=5=6 ;=1
5=6< :\$454653egw2tgg2f4 :5462



Simultaneous measurement of O and Co, Pt. with 2 beam.

CoOOH

electron hole kinetics.



Necessary conditions

Two beams have the same size to monitor the same position,

Two beams have the same timing after the pump beam.

The same pulse width